

FireSilver-Staining-Protocol

Materials

Kit Contents (Stock Solutions, store at 4°C)

20 mg Sensitizer 1 (10 x)
1 x Sensitizer 2 (10 ml)
100 x Silver Solution (10 ml)
5 x Pre-Developer Solution (200 ml)
5 x Developer 1 Solution (200 ml)
100 x Developer 2 Solution (10 ml)
5 x Stop Solution (200 ml)
Staining Protocol

Additional Required Materials

96% Ethanol
Acetic acid
Ultra pure water

Prepare solutions A-G (store at 4°C, prepare fresh for best results)

- A Fixing solution:** 50% Ethanol, 10% Acetic Acid in pure water (25 ml Ethanol, 5 ml Acetic Acid adjust to 50 ml with water).
- B Ethanol Wash solution:** 20% Ethanol (10 ml Ethanol adjust to 50 ml with water).
- C Sensitizer solution:** Add 1ml **Sensitizer 2** to 20mg **Sensitizer 1** (use solution not longer than one week!), dilute this solution before use 1:100 with ultra pure water, (500µl - adjust to 50 ml with water).
- D Silver solution:** Dilute the **Silver Stock Solution** 1:100 with ultra pure water before use (500µl - adjust to 50 ml with water).
- E Pre-developer solution:** Dilute the **Pre-Developer Stock Solution** 1:5 with ultra pure water before use (10 ml - adjust to 50 ml with water).
- F Developer solution:** Dilute the **Developer 1 Stock Solution** 1:5 with ultra pure water and pipette 1/100 volume **Developer 2 Stock Solution** to the mixture (10 ml Developer stock solution 1 - adjust to 50 ml with water plus 500µl Developer Stock Solution 2).
- G Stop solution:** Dilute the **Stop Solution** 1:5 with ultra pure water before use (10 ml - adjust to 50 ml with water).

Staining procedure

1. Fixing - After electrophoresis place the gel (7x8 cm) into a clean tray with 50 ml of the Fixing solution **A** at least for 30 minutes up to over night.
Note: A Longer fixing period increase the sensitivity.
2. Ethanol wash - Decant the Fixing solution and wash the gel for 10 minutes with 50 ml of the 20% Ethanol solution **B**.
3. Sensitization – Incubate the gel for 1 minute with 50 ml of the Sensitizer solution **C**.
4. Water wash 1 – Wash the gel for 1 minute with 100 ml of ultra pure water.
5. Water wash 2 – Repeat step 4.
6. Silver equilibration – Equilibrate the gel for 30 minutes with 50 ml of the Silver solution **D**.
7. Water wash 3 – Wash the gel for 30 seconds with 50 ml of ultra pure water.
Note: Washing longer than 1.5 minutes will result in decreased sensitivity.
8. Gel pre-development – Pre-develop the gel for 1 minute with 50 ml of the Pre-developer solution **E**.
Note: Incubation longer than 2 minutes may results in increased background.
9. Gel development – Develop the gel for 1-15 minutes with 50 ml of the Developer solution **F***.
Note: Development times longer as 15 minutes may be required to detect bands or spots with very low protein concentrations ($< 1 \text{ ng/mm}^2$).
10. Stop – Incubate the gel for 20 minutes in Stop solution **G**.
11. Storage – Store the gel in ultra pure water.

Note: Shake the gel during the process gently.
